

Determination of Deoxynivalenol, Zearalenone, T-2/HT-2 toxins, Fumonisin B1 and B2 in maize by multiplex dipstick immunoassay

MycoRed Training Course

Dr. N. Nivarlet

October 6th 2010



UNISENSOR

Food Diagnostic
Engineering

unisensor 
www.unisensor.be

Supplying Innovative and Advanced Products for
the Rapid Analysis of Food Products.

Organisation

- From R&D to Market
- R&D partner in EU-FP7 (CONFIDENCE)
- International distribution network
- A team of 23 people

Products

- Antibiotics (TetraSensor, TwinSensor, SulfaSensor, QuinoSensor, CapSensor)
- Authenticity : CowSensor, GoatSensor, ChickenSensor, PorkSensor, BeefSensor
- Adulteration : MelamineSensor
- Mycotoxins : MycoSensor, AflaSensor


trisensor


Tetra sensor


Twin sensor ^{BT}


myco sensor


afla sensor


mela sensor fluo

Plan

- General introduction of the *Confidence* project
- Our Objectives in the *Confidence* project
- Introduction on dipstick tests
- Description of the dipstick test for **Maize** (in Lab)
 - Protocol
 - Results interpretation
- Future Plans



Plan

- General introduction of the *Confidence* project
- Our Objectives in the *Confidence* project
- Introduction on dipstick tests
- Description of the dipstick test for **Maize** (in Lab)
 - Protocol
 - Results interpretation
- Future Plans



Confidence

Confidence

Contaminants in food and feed : Inexpensive
Detection for Control of Exposure

- FP7 collaborative Project
- Duration : 48 months (2008-2012)
- Participants : 17 partners from 10 countries, representing universities, research institutes, industry and SMEs



The target contaminants

➤ Cluster 1 : *Organic Pollutants*

- POPs (WP1a)
- Perfluorinated compounds (PFCs) (WP1b)
- Pesticides (WP1c)

➤ Cluster 2 : *Veterinary pharmaceuticals*

- Coccidiostats (WP2a)
- Antibiotics (WP2b)

➤ Cluster 3 : *Heavy metals speciation* (WP3)

➤ Cluster 4 : *Biotoxins*

- Alkaloids (WP4a)
- Marine Biotoxins (WP4b)
- Mycotoxins (WP4c)



The commodities in *Confidence*

- Fish/shellfish, fish feed



- Cereals, cereal-based feed and food



- Potatoes/vegetables



- Honey, dairy products



- Eggs, meat



Plan

- General introduction of the *Confidence* project
- Our Objectives in the *Confidence* project
- Introduction on dipstick tests
- Description of the dipstick test for **Maize** (in Lab)
 - Protocol
 - Results interpretation
- Future Plans



Objective

Develop Multiplex dipstick tests

for the determination of the *Fusarium* toxins DON, ZEA and T-2/HT-2 in **wheat, barley and oats**;

for DON, ZEA, T-2/HT-2 and FBs in **Maize** and maize by-products for **feed**;

for DON, ZEA, T-2/HT-2 and FBs in cereal based **food**.

- By developing the multiplex dipstick in buffer
- By adapting the multiplex dipstick in real matrices




Plan

- General introduction of the *Confidence* project
- Our Objectives in the *Confidence* project
- Introduction on dipstick tests
- Description of the dipstick test for **Maize** (in Lab)
 - Protocol
 - Results interpretation
- Future Plans



Dipstick Tests

What is a dipstick?



Need of reagents :

- competitors (ISPA, Unisensor)
- immunogens (Unisensor)
- Antibodies (CER)




For each test line
(ZEA, T2, DON, FB1)



Dipstick Tests

- *Design of Multiplex dipsticks (ZEA, T2/HT2, DON, FBs)*





MULTIPLEX DIPSTICK





Dipstick Tests

- *Design of Multiplex dipsticks (ZEA, T2/HT2, DON, FBs)*





Dipstick Tests

- *Design of Multiplex dipsticks (ZEA, T2/HT2, DON, FBs)*




Dipstick Tests

- *Design of Multiplex dipsticks (ZEA, T2/HT2, DON, FBs)*



Dipstick Tests

- *Design of Multiplex dipsticks (ZEA, T2/HT2, DON, FBs)*




MULTIPLEX DIPSTICK

POSITIVE sample




Dipstick Tests

- *Design of Multiplex dipsticks (ZEA, T2/HT2, DON, FBs)*



MULTIPLEX DIPSTICK



Test Lines equal to CTRL line

LOW POSITIVE sample



Dipstick Tests


The CONTROL line will be used to validate the dipstick test and also to compare all the test lines to an internal reference.

No CTRL line = Invalid
CTRL line = Valid

In order to conclude if the test is POSITIVE or NEGATIVE, the reader will compare each test line to the CTRL line.
Since each test line doesn't give positive result for the same ratio (usually 1), the reading with the readsensor is a must



ReadSensor



Plan

- General introduction of the *Confidence* project
- Our Objectives in the *Confidence* project
- Introduction on dipstick tests
- Description of the dipstick test for **Maize** (in Lab)
 - Protocol
 - Results interpretation
- Future Plans



Multiplex dipstick tests in MAIZE

10g of ground sample



+ 40ml of Water



Blending 2 minutes

+ 60ml of Methanol



Blending 2 minutes

Dilution 10 times in Buffer 1



Run the dipstick test

E
X
T
R
A
C
T
I
O
N


Total time
4min.



Multiplex dipstick tests in MAIZE

PROTOCOL

- 1) Dilute the extract 10 times in Buffer 1
- 2) Add 200 μ l of the diluted extract in the well, mix 20 times to homogenize
- 3) Let Incubate 10 minutes @ 40°C
- 4) Add the strips in the wells and wait for another 10 minutes
- 5) Remove the filter and read the strips with the Readsensor



CUT OFF

	POS	LPOS	NEG
<u>ZEA</u>	0	1.0	1.3
<u>T2</u>	0	1.3	1.6
<u>DON</u>	0	1.0	1.3
<u>FB1</u>	0	0.3	0.5



Limits of Detection

➤ 80% of EU levels

Desired LOD (µg/kg)	DON	ZEA	T2+HT-2	FB1+FB2
Maize	1400	280	400	3200
Durum Wheat/Oats	1400	80	400	-
Maize feed (corn gluten)	9600	2400	400	4000
Maize food	400	80	80	640
Oat/Wheat food	400	40	80	-



Plan

- General introduction of the *Confidence* project
- Our Objectives in the *Confidence* project
- Introduction on dipstick tests
- Description of the dipstick test for **Maize** (in Lab)
 - Protocol
 - Results interpretation
- Future Plans



Future Plans

- Validation : Scheduled for **JULY 2011** (inside the confidence)
- Launch of the Product : “MYCOSENSOR”
Scheduled for the **beginning of 2011**
- Price : to be determined
- Distributors : Distributors of UNISENSOR’s Products
- General Informations : info@unisensor.be
- Commercial Purposes : contact olivier.heyne@unisensor.be



....That's time for the
demonstration in Lab....



Thank you !

